

**OLIVE VIEW-UCLA MEDICAL CENTER  
DEPARTMENT OF PATHOLOGY  
POLICY & PROCEDURE**

**NUMBER: 11906  
VERSION: 1**

**SUBJECT/TITLE: SS 013 - OIL RED O STAINING PROCEDURE**

**POLICY:** The Oil Red O staining procedure shall be performed at the request of the pathologist.

**PURPOSE:** To demonstrate neutral lipids in frozen tissue sections.

**DEPARTMENTS:** **PATHOLOGY & LABORATORY SERVICES**

**DEFINITIONS:**

**PROCEDURE:** **FIXATIVE**  
10% neutral buffered formalin. No alcoholic fixatives should be used.

**EQUIPMENT**

Cryostat                      Coplin jars                      Filter paper  
Erlenmeyer flasks      Graduated Cylinders  
Portable Scale w/ Weigh Boats  
Portable Stirring Plate w/ Stir Bars

**SLIDE PREPARATIONS**

Cut frozen sections at 10 microns. Paraffin sections **cannot** be utilized as dehydrating and clearing agents dissolve the fat.

**QUALITY CONTROL**

Most tissue contains some fat, but a touch prep made with fatty tissue can also be used.

**REAGENTS**

**I. Stock Solutions**

- Oil Red O Stock Solution
  - Oil Red O.....2.5 gm
  - 99% Isopropanol Alcohol.....500 ml

Mix Well

- 60% Isopropanol Alcohol
  - Isopropanol Alcohol.....60 ml
  - Distilled Water.....40 ml

**SUBJECT/TITLE: SS 013 - OIL RED O STAINING PROCEDURE**

**Policy Number: 11906**

**Page Number: 2**

**II. Working Solutions**

- Oil Red O Stock Solution.....24 ml
- Distilled Water.....16 ml

Mix well and allow to stand for 10 minutes and afterwards filter the solution.

**III. Pre-made Solutions**

- Lithium Carbonate
- Gill Hematoxylin

**PROCEDURE**

1. Cut frozen sections at 10 microns. Fix in 37-40% formaldehyde (conc.) for 1 minute.
2. Wash well in tap water. Blot off excess water.
3. Stain sections in oil red O for 10 minutes.
4. Wash sections in tap water. Check under the microscope.
5. Differentiate briefly in 60% Isopropanol if fat is smeared on top of section.
6. Wash sections in tap water.
7. Counterstain sections in Gill hematoxylin for 1 minute.
8. Wash sections well in tap water.
9. Blue in Lithium Carbonate.
10. Wash sections well in tap water.
11. Mount section with an aqueous mounting medium.

**Note:** The fat in the section is relatively liquid and mobile, so care should be taken that no pressure is placed on the coverglass or the fat may be displaced. If air bubbles are present in the section, remove the coverglass by soaking the slide in warm water.

**RESULTS:**

Fat..... Intense Red  
Nuclei .....Blue

References: Freida L. Carson, PhD, HT (ASCP), Histotechnology, A Self-Instructional Text, ASCP Press, ASCP , Chicago, 1997	
Approved by: Holli Mason (Laboratory Director)	Date: 05/22/2020
Review Dates: 05/22/2020	
Next Review: 05/22/2022	
Distribution: Pathology & Laboratory Services	
Original Date: 10/1/1992	

**SUBJECT/TITLE: SS 013 - OIL RED O STAINING PROCEDURE**

**Policy Number: 11906**

**Page Number: 3**