



TO PRINT A MICROBIOLOGY WORKLIST IN CERNER




1. Open up Microbiology Result Entry




2. Click on Workqueue.
3. Click on Create Workqueue.



4. Under Report, choose Procedure Worklist. Then, scroll down to the HDH choices under Display. Choose one of the HDH Display Sequences. Click on it to highlight it, then click on OK.



5. The Create Workqueue box pops up.
6. Choose Create Queue.



7. Once the Queue is created, the Report Output – preview box will appear. Click on the print icon.
8. Choose the correct printer, and click OK.
9. Once printed, close this Workqueue box, so a new queue can be created.

