

**Warming Cabinet for Solutions and Blankets
ATTACHMENT 1
Solution/Linen Warming Cabinet Temperature Log**

INSTRUCTIONS:

- A. Record the temperature of the warmer daily in the appropriate section next to the date.
Note: Acceptable maximum temperature: 40°C for IV Solutions and 54°C F for Linens
- B. Adjust temperature as needed to maintain acceptable range.
- C. Indicate corrective actions or changes to adjust temperature in the comments section of log.
Note: Report malfunction of warmer to Facilities Management Department x73991
- D. Begin a new form on the first of each month and return completed logs to Nurse Manager.

Month: ____ Year: ____ Location of Warmer: _____ Type of Warmer

(Top): IV solution Irrigation solution (Bottom): Linen

Date	Temperature		Initial: temp check and expired IV bag check completed	Comments/Actions
	Top	Bottom		
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				

If malfunctioning call facilities #73991 – Discard fluids if temp exceeds 40°C and if fluids are not used within the 14-day warmed shelf life.